



Surface Des Plus

VAH listed

Concentrated disinfectant and cleaner for surfaces

- free from aldehyde
- strong cleaning properties
- suitable for PVC and linoleum floorings



Product properties

Adelhyde-free concentrate for disinfecting and cleaning of medical products, medical inventory and surfaces of all sorts.

Composition

100 g solution contain:
16 g Alkyl(C12-16)dimethylbenzyl ammonium chloride.

Chemical physical data

Appearance of concentrate: clear, blue-green solution
pH-value: 8,5 +/-0,5
Density: approx. 1,00 g/ml

Range of Application

Acc. to Medical Devices Act:
Aldehyde-free disinfection and cleaning of medical devices such as respiratory masks as well as incubators.

Acc. to Biocidal Products Regulation:
Disinfection and cleaning of all kinds of surfaces and medical inventory.

Direction for use

Surface disinfection by wiping: Apply a sufficient quantity of the prepared Safeline Surface Des plus working solution to the surfaces to be disinfected. Use mechanical action to completely wet all surfaces (scrub-wipe disinfection).

Concentration and exposure times

Bactericidal incl. MRSA, yeasticidal acc. to VAH-guidelines	1% (10g/l)	5 Min.
	0.75% (7.5g/l)	15 Min.
	0.5% (5g/l)	60 Min.
tuberculocidal	0.5% (5g/l)	60 Min.
HBV/HIV/HCV/ BVDV, Vaccinia (acc. to the guidelines of the Robert Koch Institute)	0.5% (5g/l)	15 Min.
Noro virus (MNV)	1% (10g/l)	60 Min.
Papova/Polyomaviren/SV40	1% (10g/l)	60 Min.

Delivery Forms

Single package	Sales unit	Ref.-No.
5 Litre canister	1 x 5 L	09950

Pharmaceutical expert report

Prof. Dr. Schubert, Frankfurt/Main:
Certificate for Surface Disinfection according to DGHM/VAH Guidelines from 15.04.2010 Certificate for Tbc effectiveness from 15.04.2010

Dr. med. F.-A. Pitten, Hygieniker, Gießen:
Certificate for Surface Disinfection according to DGHM/VAH Guidelines from 14.04.2010

Prof. Dr. med. H.-P. Werner, Bischofshofen:
Certificate for Surface Disinfection according to DGHM/VAH Guidelines from 28.04.2014

Classification according regulation (EC) Nr. 1272/2008:



Danger:
n-alkylbenzyl-di-methyl ammonium chloride; alkylpropylendiamine-1,5-bis-guanidinium acetate; tridecylamine.
Causes severe skin burns and eye damage. Very toxic to aquatic life with long lasting effects. Wear protective gloves/ protective clothing/eye protection/ face protection. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/ shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Call a POISON CENTER/doctor. Dispose of contents/container to approved disposal company or local collection.

Use disinfectants in a safe way. Always read label and product information before use.

AMPri

AMPri Handelsgesellschaft mbH
Benzstr. 16 / 21423 Winsen (Luhe), Germany
Tel. +49 (0)4171 - 8480-0 / www.ampri.de / info@ampri.de

bioSight GmbH
Brasselsbergstr. 3 / 34132 Kassel, Germany
Tel.: +49 (0)561 - 81678111