

### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

# **Technical Data Sheet**

Article-No.: 05012-W

EN

Description: **MED-COMFORT** 

non-woven visitor-coat

with stand-up-collar and zip



SIZES				
	M	L	XL	XXL
width	63 cm	66 cm	69 cm	73 cm
length	83 cm	85 cm	90 cm	92 cm

PRODUCT DES	CR	IPTION										
material	•	Polypro- pylene (PP)		Polyethy- lene (PE)		SMS						
colour		blue [B]		green [G]		orange [O]	red [R]	5	white [W]	yellow [Y]		light-blue [LB]
		dark-blue [DB]		dark-green [DG]								
closure		to tie	Z	zipper		push bottons	velcro					
characteristics		tricot-cuff	Į.	elastic wrist- cuff	V	latex-free	hood	v	stand-up collar	partially laminated		sewed PE- oversleeves
		shortenable by perforation		blocked		dust-cover	embossed		short sleeves	with cardboard	v	single folded

<b>REGULATORY</b>	AF	FAIRS												
PPE-Regulation (EU) 2016/425		category I		category II		category III	V	no PPE- article						
MD-Regulation (EU) 2017/745		class I		class II		class III		sterile		measuring fundtion	7	no medical device		
Food Contact (EG) 1935/2004	7	acidic foods	>	aqueous foods	>	fatty foods	V	alcoholic foods	>	dry foods		not approved for food- contact	5	ור. ור

STANDARDS							
tested in	□ DIN EN	□ ISO 16603	□ ISO 16604	□ ISO 22610	□ ISO 22611	□ ISO 22612	
accordance to	14126						
	□ EN ISO 13688	□ EN 13034					

1/3 Herausgabe: 15.09.2021



### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

## **Technical Data Sheet**

Article-No.: 05012-W

Description: MED-COMFORT

non-woven visitor-coat

with stand-up-collar and zip



LOGISTIC DATA SUBPACKING					
general information					
material		bag			
pieces per subpacking	5	1			
GTIN subpacking size	М	4044941020594			
GTIN subpacking size	L	4044941020563			
GTIN subpacking size	XL	4044941020570			
GTIN subpacking size	XXL	4044941020587			
measures & sizes					
length		270 mm			
width		340 mm			
height		-			
weights					
size	net weight	brut weight			
M	71 g	76 g			
L	73 g	78 g			
XL	77 g	82 g			
XXL	80 g	85 g			

LOGISTIC DATA PA				
general information				
kind of palette		euro-palette		
measures & sizes				
cartons per layer		5		
layer per palette		4		
palette height		119		
weights				
size	net weight	brut weight		
M	162 g	187 g		
L	166 g	191 g		
XL	174 g	199 g		
XXL	180 g	205 g		

EN

LOGISTIC DATA OUTER PACKING					
general information					
material		carton			
subpackings per outer	packing	100			
GTIN outer packing siz	e M	4044941020594			
GTIN outer packing siz	e L	4044941020600			
GTIN outer packing siz	e XL	4044941020617			
GTIN outer packing siz	e XXL	4044941020624			
measures & sizes					
length		500 mm			
width		300 mm			
height		260 mm			
weights					
size	net weight	brut weight			
М	7.600 g	8.100 g			
L	7.800 g	8.300 g			
XL	8.200 g	8.700 g			
XXL	8.500 g	9.000 g			



### AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

## **Technical Data Sheet**

Article-No.: 05012-W

Description: **MED-COMFORT** 

> non-woven visitor-coat with stand-up-collar and zip



### WARNINGS AND SAFETY INFORMATION

storage / expiry	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms,
	as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The
	expiry date - valid for proper storage - is stated on the packaging.

use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product	
	is suitable for the intended use.	
	Ensure that the size matches that of the user. Do not make any changes to the product.	
	Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).	
	This product is intended for single use only	

disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company.
	Unused products can be disposed of with household waste.

warnings	Prolonged wearing can lead to heat stress.	
	This product is flammable - keep away from fire.	

#### donning and doffing WAY OF PUTTING ON:

instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body

body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 6

date: 21. 10.2024 changes and errors excepted

QMFORM\_67.001 Herausgabe: 15.09.2021