

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05070-W-BLOCK

05070-B-BLOCK

MED-COMFORT Description:

Uni

SIZES

disposable PE apron

medium, smooth, blocked, in a bag



width	/5 cm						
length	120 cm						
PRODUCT DES	CRIPTION						
material	Polypro-	Polyethy- lene (PE)	☐ SMS				
colour	✓ blue [B]	☐ green [G]	orange [O]	☐ red [R]	white [W]	☐ yellow [Y]	☐ light-blue [LB]
	☐ dark-blue [DB]	☐ dark-green [DG]					
closure	✓ to tie	☐ zipper	push bottons	☐ velcro			
characteristics	☐ tricot-cuff	elastic wrist- cuff	☑ latex-free	hood	stand-up collar	partiallylaminated	sewed PE- oversleeves
	shortenable by perforation	☑ blocked	□ dust-cover	☐ embossed	☐ short sleeves	☐ with cardboard	□ single folded
	•		•		•	•	•
REGULATORY	AFFAIRS						
PPE-Regulation		category II	category III	☐ no PPE-			

REGULATORY A	AFFAIRS						
PPE-Regulation		category II	category III	no PPE-			
(EU) 2016/425				article			
							CE
MD-Regulation	class I	☐ class II	class III	☐ sterile	☐ measuring	✓ no medical	
(EU) 2017/745					fundtion	device	
Food Contact	□ acidic foods	☑ aqueous		☑ alcoholic	☑ dry foods	☐ not	
(EG) 1935/2004		foods		foods		approved	רא
						for food-	
						contact	
							<u> </u>

STANDARDS							
tested in	☐ DIN EN	☐ ISO 16603	☐ ISO 16604	☐ ISO 22610	☐ ISO 22611	☐ ISO 22612	
accordance to	14126						
	■ EN ISO	☐ EN 13034					
	13688						

QMFORM_67.001 EN 1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: **05070-W-BLOCK**

05070-B-BLOCK

Description: MED-COMFORT

disposable PE apron

medium, smooth, blocked, in a bag



LOGISTIC DATA SUBPACKING				
general information				
material		bag		
pieces per subpacking	5	50		
GTIN subpacking blue	· [B]	4044941020693		
GTIN subpacking whit	e [W]	4044941020105		
measures & sizes				
length		480 mm		
width		320 mm		
height		-		
weights				
size	net weight	brut weight		
uni	1.590 g	1.595 g		

LOGISTIC DATA PALETTE				
general information				
kind of palette		euro-palette		
measures & sizes				
cartons per layer 6				
layer per palette 6				
palette height	150			
weights				
size	net weight	brut weight		
uni 592 g		617 g		

EN

LOGISTIC DATA OUTER PACKING				
general information				
material		carton		
subpackings per outer	packing	10		
GTIN outer packing blo	ue [B]	4044941020709		
GTIN outer packing wl	4044941020280			
measures & sizes				
length 390 mm				
width	300 mm			
height	225 mm			
weights				
size	net weight	brut weight		
uni	16.450 g			

2/3



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05070-W-BLOCK

05070-B-BLOCK

MED-COMFORT Description:

disposable PE apron

medium, smooth, blocked, in a bag



WARNINGS AND SAFETY INFORMATION

storage / expiry	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do
date	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms,
	as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The
	expiry date - valid for proper storage - is stated on the packaging.

use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use.
	Ensure that the size matches that of the user. Do not make any changes to the product.
	Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).
	This product is intended for single use only

disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company.	l
	Unused products can be disposed of with household waste.	

warnings	Prolonged wearing can lead to heat stress.
	This product is flammable - keep away from fire.

donning and doffing

instructions

WAY OF PUTTING ON:

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 6

22. 10.2024 date:

changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021