

Technical Data Sheet

Article-No.: 05015

Description:

MED-COMFORT non woven protective coat with sewed-on apron and oversleeves made of PE



SIZES	
	Uni
width	150 cm
length	150 cm

PRODUCT DES	CR	IPTION								
material	2	Polypro- pylene (PP)	Polyethy- lene (PE)		SMS					
colour	•	blue [B]	green [G]		orange [O]	red [R]	white [W]	yellow [Y]		light-blue [LB]
		dark-blue [DB]	dark-green [DG]							
closure	7	to tie	zipper		push bottons	velcro				
characteristics		tricot-cuff	elastic wrist- cuff	-	latex-free	hood	stand-up collar	partially laminated	2	sewed PE- oversleeves
		shortenable by perforation	blocked		dust-cover	embossed	short sleeves	with cardboard	7	single folded

REGULATORY	AFFAIRS						
PPE-Regulation (EU) 2016/425	category I	category II	category III	no PPE- article			(6
MD-Regulation (EU) 2017/745	🗆 class I	class II	□ class III	□ sterile	measuring fundtion	 no medical device 	
Food Contact (EG) 1935/2004	☑ acidic foods	 aqueous foods 	fatty foods	alcoholic foods	☑ dry foods	 not approved for food- contact 	Ĩ

STANDARDS							
tested in	DIN EN	ISO 16603	ISO 16604	ISO 22610	ISO 22611	ISO 22612	
accordance to	14126						
	EN ISO	EN 13034					
	13688						



AMPri Handelsgesellschaft mbH Benzstraße 16 21423 Winsen (Luhe) - Germany Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05015

Description:

MED-COMFORT non woven protective coat with sewed-on apron and oversleeves made of PE



LOGISTIC DATA SUBPACKING							
general information							
material		bag					
pieces per subpacking	5	10					
GTIN subpacking		4044941020983					
measures & sizes							
length		480 mm					
width		320 mm					
height		-					
weights	weights						
size	brut weight						
150x150 cm 690 g 700 g							

LOGISTIC DATA PALETTE						
general information						
kind of palette euro-palette						
measures & sizes						
cartons per layer 4						
layer per palette 6						
palette height 207						
weights						
size	net weight	brut weight				
150x150 cm 180 g 205 g						

LOGISTIC DATA OUTER PACKING

general information						
material		carton				
subpackings per outer	[,] packing	10				
GTIN outer packing		4044941020990				
measures & sizes						
length		560 mm				
width		280 mm				
height		320 mm				
weights						
size	net weight	brut weight				
150x150 cm 7.000 g 7.500 g						



Technical Data Sheet

Article-No.: 05015

Description:

MED-COMFORT non woven protective coat with sewed-on apron and oversleeves made of PE



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. Ensure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). This product is intended for single use only
disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
warnings	Prolonged wearing can lead to heat stress. This product is flammable - keep away from fire.
donning and doffin instructions	 WAY OF PUTTING ON: Make sure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body body and fastening it with a bow. REMOVAL METHOD: Undo the waist and collar ties and remove the gown.

changes and errors excepted