

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05060-B-Block / 05060-W-Block

Description: **MED-COMFORT**

Uni

75 cm

SIZES

width

EN

disposable PE apron

medium, smooth, blocked, in a bag



leligui	120 (111						
PRODUCT DES	CRIPTION						
material	Polypro-pylene (PP)	Polyethy- lene (PE)	☐ SMS				
colour	✓ blue [B]	green [G]	☐ orange [O]	☐ red [R]	white [W]	☐ yellow [Y]	☐ light-blue [LB]
	☐ dark-blue [DB]	☐ dark-green [DG]					
closure	✓ to tie	☐ zipper	☐ push bottons	□ velcro			
characteristics	☐ tricot-cuff	elastic wrist- cuff	✓ latex-free	hood	stand-up collar	partiallylaminated	sewed PE- oversleeves
	shortenable by perforation	☑ blocked	☐ dust-cover	☐ embossed	□ short sleeves	with cardboard	☐ single folded
2501114502	A =						
REGULATORY PPE-Regulation	AFFAIRS category I	☐ category II	☐ category III	□ no PPE-			
(EU) 2016/425	Category	category ii	category iii	article			C€
MD-Regulation (EU) 2017/745	☐ class I	☐ class II	☐ class III	☐ sterile	☐ measuring fundtion	ono medical device	
Food Contact (EG) 1935/2004	∡ acidic foods	✓ aqueous foods	☑ fatty foods	✓ alcoholic foods	dry foods	not approved for food- contact	7.7
STANDARDS tested in	DIN EN	□ ISO 16603	ISO 16604	□ ISO 23610	ISO 22611	☐ ISO 22612	
accordance to	14126	☐ ISO 16603	☐ ISO 16604	☐ ISO 22610	☐ ISO 22611	130 22612	
	☐ EN ISO 13688	☐ EN 13034					

QMFORM_67.001 1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05060-B-Block / 05060-W-Block

Description: MED-COMFORT

disposable PE apron

medium, smooth, blocked, in a bag



LOGISTIC DATA SUBPACKING			
general information			
material	bag		
pieces per subpacking	50		
GTIN subpacking blue	4044941020082		
GTIN subpacking whit	4044941020099		
measures & sizes			
length	480 mm		
width	320 mm		
height	-		
weights			
size	net weight	brut weight	
uni	815 g	820 g	

LOGISTIC DATA PALETTE				
general information				
kind of palette		euro-palette		
measures & sizes				
cartons per layer	6			
layer per palette		10		
palette height	150			
weights				
size	net weight	brut weight		
uni	522 g	547 g		

LOGISTIC DATA OUTER PACKING			
general information			
material	carton		
subpackings per outer	10		
GTIN outer packing blu	4044941020273		
GTIN outer packing wh	4044941020266		
measures & sizes			
length	390 mm		
width	300 mm		
height	135 mm		
weights			
size	net weight	brut weight	
uni	8.200 g	8.700 g	

QMFORM_67.001 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05060-B-Block / 05060-W-Block

Description: MED-COMFORT

disposable PE apron

medium, smooth, blocked, in a bag



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. Ensure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). This product is intended for single use only

disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company.
	Unused products can be disposed of with household waste.

warnings	Prolonged wearing can lead to heat stress.
	This product is flammable - keep away from fire.

donning and doffing WAY OF PUTTING ON:

instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 5

date: 22. 10.2024

changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021