

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Uni

100 cm

SIZES

width

Article-No.: **05071-B-BLOCK-150**Description: **MED-COMFORT**

disposable PE apron

medium, smooth, blocked, in a bag



| WIGHT | 100 cm | | | | | | |
|--------------------------------|----------------------------|------------------------|----------------|-------------------|-----------------|--------------------------------|--------------------------|
| length | 150 cm | | | | | | |
| | | | | | | | |
| PRODUCT DE | SCRIPTION | | | | | | |
| material | ☐ Polypro- | ✓ Polyethy- | ☐ SMS | | | | |
| | pylene (PP) | lene (PE) | | | | | |
| colour | ✓ blue [B] | green [G] | orange [O] | ☐ red [R] | ☐ white [W] | yellow [Y] | ☐ light-blue [LB] |
| | ☐ dark-blue | ☐ dark-green | | | | | |
| | [DB] | [DG] | | | | | |
| closure | ✓ to tie | ☐ zipper | push bottons | □ velcro | | | |
| characteristics | ☐ tricot-cuff | elastic wrist- cuff | ☑ latex-free | hood | stand-up collar | ☐ partially laminated | sewed PE- oversleeves |
| | shortenable by perforation | ☑ blocked | □ dust-cover | ☐ embossed | ☐ short sleeves | ☐ with cardboard | ☐ single folded |
| REGULATORY | AFFAIRS | | | | | | |
| PPE-Regulation | | ☐ category II | ☐ category III | no PPE- | | | |
| (EU) 2016/425 | | | | article | | | CE |
| MD-Regulation | ☐ class I | class II | class III | sterile | ☐ measuring | ✓ no medical | 66 |
| (EU) 2017/745 | | | | | fundtion | device | |
| Food Contact (EG) 1935/2004 | ✓ acidic foods | ✓ aqueous foods | ✓ fatty foods | ✓ alcoholic foods | ☑ dry foods | not approved for food- contact | 77 |
| STANDARDS | | | | | | | |
| tested in | ☐ DIN EN | ISO 16603 | ☐ ISO 16604 | ISO 22610 | ☐ ISO 22611 | ISO 22612 | |
| accordance to | 14126 | .50 10003 | 130 10034 | .55 22510 | | | |
| | ☐ EN ISO 13688 | ☐ EN 13034 | | | | | |

QMFORM_67.001 1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: **05071-B-BLOCK-150**

Description: MED-COMFORT

disposable PE apron

medium, smooth, blocked, in a bag



| LOGISTIC DATA SUBPACKING | | | |
|--------------------------|---------------|-------------|--|
| general information | | | |
| material | bag | | |
| pieces per subpacking | 50 | | |
| GTIN subpacking | 4044941021737 | | |
| measures & sizes | | | |
| length | | 600 mm | |
| width | 320 mm | | |
| height | - | | |
| weights | | | |
| size | net weight | brut weight | |
| uni | 2.655 g | 2.660 g | |

| LOGISTIC DATA PALETTE | | | |
|------------------------------|------------|-------------|--|
| general information | | | |
| kind of palette euro-palette | | | |
| measures & sizes | | | |
| cartons per layer 5 | | | |
| layer per palette | 10 | | |
| palette height | 130 | | |
| weights | | | |
| size | net weight | brut weight | |
| uni | 557 g | 582 g | |

| LOGISTIC DATA OUTER PACKING | | | |
|-----------------------------|------------|---------------|--|
| general information | | | |
| material | | carton | |
| subpackings per outer | packing | 4 | |
| GTIN outer packing | | 4044941021744 | |
| measures & sizes | | | |
| length | | 525 mm | |
| width | | 300 mm | |
| height | | 115 mm | |
| weights | | | |
| size | net weight | brut weight | |
| uni 10.640 g 11.140 g | | 11.140 g | |



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05071-B-BLOCK-150 Description: **MED-COMFORT**

disposable PE apron

medium, smooth, blocked, in a bag



WARNINGS AND SAFETY INFORMATION

| storage / expiry | Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| date | not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging. |
| use and control | Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. |
| | Ensure that the size matches that of the user. Do not make any changes to the product. |
| | Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). |
| | This product is intended for single use only |
| | |
| diamanal | |

| disposal | The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. |
|----------|-----------------------------------------------------------------------------------------------------------------------|
| | Unused products can be disposed of with household waste. |

| warnings | Prolonged wearing can lead to heat stress. |
|----------|--------------------------------------------------|
| | This product is flammable - keep away from fire. |
| | |
| | |
| | |
| | |
| | |
| | |

donning and doffing WAY OF PUTTING ON:

instructions

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).
- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 3

date: 22. 10.2024 changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021