

AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: **05500**

Description: MED-COMFORT

non-woven visitor-coat

light, size M/L



SIZES	
	M/L
width	135 cm
length	115 cm

EN

PRODUCT DESC	CRIPTION						
material	✓ Polypro-	☐ Polyethy-	☐ SMS				
	pylene (PP)	lene (PE)					
colour	☐ blue [B]	green [G]	orange [O]	red [R]	white [W]	yellow [Y]	☐ light-blue
							[LB]
	✓ dark-blue	☑ dark-green					
	[DB]	[DG]					
closure	✓ to tie	☐ zipper	☐ push	☐ velcro			
			bottons				
characteristics	☐ tricot-cuff	elastic wrist-	☐ latex-free	☐ hood	stand-up	☐ partially	sewed PE-
		cuff			collar	laminated	oversleeves
	□ shortenable	□ blocked	☐ dust-cover	□ embossed	□ short sleeves	with	
	by	blocked	uust-covei	embossed	□ SHOLL Sleeves	cardboard	Siligle folded
	perforation					caraboara	
		I.	I.		I	I.	I.
REGULATORY A	AFFAIRS						
PPE-Regulation	☐ category I	☐ category II	☐ category III	✓ no PPE-			
(EU) 2016/425				article			
MD-Regulation	☐ class I	class II	☐ class III	sterile	☐ measuring		
(EU) 2017/745					fundtion	device	
Food Contact	acidic foods	☑ aqueous		☐ alcoholic	☑ dry foods	□ not	
(EG) 1935/2004	delate loods	foods	latty 10003	foods	ary roous	approved	[]"
(==, ====, === :						for food-	K
						contact	
				•	•		•
STANDARDS							
tested in	☐ DIN EN	☐ ISO 16603	☐ ISO 16604	☐ ISO 22610	☐ ISO 22611	☐ ISO 22612	
accordance to	14126						
	- FNUSO	□ FN 4303 *					
	☐ EN ISO 13688	☐ EN 13034					
	13000						

QMFORM_67.001 1/3 Herausgabe: 15.09.2021



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de


Technical Data Sheet

Article-No.: 05500

Description: MED-COMFORT

non-woven visitor-coat

light, size M/L



LOGISTIC DATA PALETTE			
general information			
kind of palette	kind of palette euro-palette		
measures & sizes			
cartons per layer		7	
layer per palette		7	
palette height	190		
weights			
size	net weight	brut weight	
M/L	159 g	184 g	



LOGISTIC DATA OU	TER PACKING		
general information			
material		carton	
subpackings per outer	packing	10	
GTIN outer packing blu	ue [B]	-	
GTIN outer packing gre	een [G]	-	
GTIN outer packing or	ange [O]	-	
GTIN outer packing re-	-		
GTIN outer packing wh	nite [W]	4044941021416	
GTIN outer packing ye	llow [Y]	-	
GTIN outer packing lig	ht-blue [LB]	-	
GTIN outer packing da	4044941021393		
GTIN outer packing dark-green [DG]		4044941021409	
measures & sizes			
length	460 mm		
width	220 mm		
height		250 mm	
weights			
size	net weight	brut weight	
M/L	2.750 g	3.250 g	



AMPri Handelsgesellschaft mbH

Benzstraße 16

21423 Winsen (Luhe) - Germany

Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: 05500

Description: MED-COMFORT

non-woven visitor-coat

light, size M/L



WARNINGS AND SAFETY INFORMATION

storage / expiry	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do		
	not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.		

use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use.
	Ensure that the size matches that of the user. Do not make any changes to the product.
	Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).
	This product is intended for single use only

disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.		
warnings	Prolonged wearing can lead to heat stress.		
	This product is flammable - keep away from fire.		

This product is flammable - keep away from fire.

instructions -

donning and doffing | WAY OF PUTTING ON:

- Make sure that the size matches that of the user. Do not make any changes to the product.
- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).

3/3

- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband around your body and fastening it with a bow.

REMOVAL METHOD:

- Undo the waist and collar ties and remove the gown.

rev-no.: 5

date: 22. 10.2024

changes and errors excepted

QMFORM_67.001 Herausgabe: 15.09.2021