

Technical Data Sheet

Article-No.: **02501**

Description: **HAPPY KIDS**
facemask

3-ply, with elastic-band



SIZES	
length	width
175 mm	95 mm

PRODUCT DESCRIPTION							
material	<input checked="" type="checkbox"/> non-woven PP	<input type="checkbox"/> paper					
colour	<input type="checkbox"/> blue [B]	<input type="checkbox"/> berry [Berry]	<input type="checkbox"/> bordeaux [BR]	<input type="checkbox"/> cedro [C]	<input type="checkbox"/> green [G]	<input type="checkbox"/> lagoon [Lagoon]	<input type="checkbox"/> purple [L]
	<input type="checkbox"/> lime [Lime]	<input type="checkbox"/> Mix (blue,yellow,	<input type="checkbox"/> orange [O]	<input type="checkbox"/> pink [P]	<input type="checkbox"/> white [W]	<input type="checkbox"/> yellow [Y]	<input type="checkbox"/> platinum / grey [PT]
	<input type="checkbox"/> blue patterned [DB]	<input type="checkbox"/> wine-red patterned [WR]	<input type="checkbox"/> black [SCH]	<input type="checkbox"/> black-red-gold	<input type="checkbox"/> clean ocean / ocean green [CO]	<input checked="" type="checkbox"/> White with bear print	
characteristics	<input checked="" type="checkbox"/> type II	<input checked="" type="checkbox"/> 3 layers	<input checked="" type="checkbox"/> latex-free	<input type="checkbox"/> with valve	<input type="checkbox"/> anti-fog	<input checked="" type="checkbox"/> with ear loops	<input type="checkbox"/> with activated carbon filter
characteristics	bacterial filtration capacity \geq 98%						

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	<input type="checkbox"/> category I	<input type="checkbox"/> category II	<input type="checkbox"/> category III	<input checked="" type="checkbox"/> no PPE-article			CE
MD-Regulation (EU) 2017/745	<input checked="" type="checkbox"/> class I	<input type="checkbox"/> class II	<input type="checkbox"/> class III	<input type="checkbox"/> sterile	<input type="checkbox"/> measuring fundtion	<input type="checkbox"/> no medical device	

STANDARDS							
tested in accordance to	<input checked="" type="checkbox"/> EN 14683:2019	<input type="checkbox"/> EN 149:2001					

Technical Data Sheet

Article-No.: **02501**
Description: **HAPPY KIDS**
facemask
3-ply, with elastic-band



LOGISTIC DATA SUBPACKING		
general information		
material	carton	
pieces per subpacking	50	
GTIN subpacking blue [B]	4044941031767	
measures & sizes		
length	180 mm	
width	100 mm	
height	90 mm	
weights		
size	net weight	brut weight
uni	150 g	185 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	8	
layer per palette	8	
palette height	175	
weights		
size	net weight	brut weight
uni	150 g	175 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	10	
GTIN outer packing blue [B]	4044941031774	
measures & sizes		
length	520 mm	
width	120 mm	
height	200 mm	
weights		
size	net weight	brut weight
uni	1.850 g	2.350 g

Technical Data Sheet

Article-No.: **02501**
Description: **HAPPY KIDS**
facemask
3-ply, with elastic-band



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store face mask in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Always use face masks only for their intended use. A check/risk assessment must be carried out to ensure that the face mask is suitable for the intended use. Change the mask when it is soaked through.
disposal	Used face masks must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused face masks can be disposed of with household waste.
warnings	Change the mask when it is soaked through. Article is highly flammable - keep away from heat and/or fire.
donning and doffing instructions	<p>PUTTING ON: Place the mask over the mouth and nose on the face and secure it in place with straps or ear loops. Adjust the flexible nose clip to the shape of the face so that as little air as possible can escape from the upper end.</p> <p>PULL OUT: Loosen the straps to remove the face mask from the face. Do not touch the fleece material. Disinfect your hands after disposing of the face mask.</p>