

TECHNICAL DATA SHEET



AMPri Handelsges. mbH • Benzstraße 16 • 21423 Winsen/Luhe, Germany • Tel. 04171 84 80-0 • www.ampri.de • info@ampri.de

Kidney Dish, Article No.: 7100040



Material:	stainless steel
Sterilization	Steam heat" in the autoclave Sterilize the kidney dish with moist heat, at 134°C for 4 minutes at 2 bar pressure. It is recommended that the sterilization temperatures do not exceed 137°C.
Sizes:	250 mm in length, 35 mm deep
Packing:	100 pcs per carton

07/24 Rev2